



Harvard Polycarboxylat Cement

Zinc polycarboxylate cement for permanent luting and lining

Indications

- Zinc polycarboxylate cement for permanent luting of crowns and bridges, inlays and onlays on natural core and for luting supraconstruction on implants
 - For crowns and bridges made of
 - Zirconia
 - Aluminum oxide
 - Lithiumdisilicate
 - Silicate ceramic
 - For all classical metal and metal-ceramic materials
- Liner for all types of filling materials

Properties

- Better biocompatibility than phosphate cements
- Non irritant for sensitive teeth
- Less irritant to the pulp than zinc phosphate cement
- Cement according to DIN EN ISO 9917-1

Other advantages

- Easy and safe application
- Unmatched price-performance ratio for luting materials

Article

Order no.

Order no.

Powder shade

3 - yellowish white

4 - light yellow

Clinic powder 100 g

7031203

7031204

Single powder 35 g

7031503

7031504

Liquid

Clinic liquid 40 ml

7031300

Single liquid 15 ml

7031600